


## Abstract

MuRGE (Mutualism in a Reduced Gravity Environment) is a NASA flight-research experiment to investigate the microgravity effects associated with cell-cell communication and beneficial microbe-host interactions using a plant-fungal model system. This investigation will use a clinostat, an instrument that slowly rotates the plants to negate the effects of gravitational pull on plant growth (gravitropism) and development, to simulate microgravity. I will be using the endophytic fungus *Piriformospora indica* (Pi) and the model plant species *Arabidopsis thaliana* (At). *P. indica* has been shown to colonize roots of various plant species, including *A. thaliana*, and to increase plant growth and resistance to stress. The fungus has the ability to grow from spores or in axenic cultures without the presence of a host. *P. indica* spores and *P. indica* extract will be used to inoculate *Arabidopsis* seeds germinated on a clinostat in order to determine if simulated microgravity affects the interaction between the fungus and its plant host.


















# NASA INSPIRE

Karishma Patel  
Pre-college Intern  
Kennedy Space Center, Fl  
Summer 2010




# About Me

- ▶ Karishma Patel
- ▶ Hometown: Calhoun, Georgia
- ▶ Georgia Institute of Technology, Class of 2014
- ▶ Hobbies: Running, playing the trumpet, hiking





# Overview

- ▶ Work Overview
  - ▶ MuRGE
  - ▶ Exposure to
  - ▶ Work Related Activities
  - ▶ Non-Work Related Activities
  - ▶ Lessons Learned
  - ▶ Lasting Impact
  - ▶ Future Plans
  - ▶ Acknowledgements
- 



# Work Overview

- ▶ Mentor: Michael Roberts
- ▶ Space Life Sciences Lab (SLSL)
- ▶ Dynamac Corporation
- ▶ NASA Exploration Life Support Program
- ▶ NASA-KSC Engineering Directorate (NE-S)







# MuRGE

## (Mutualism in a Reduced Gravty Environment)

- ▶ Mutualism is the way two organisms biologically interact where each individual derives a fitness benefit.
- ▶ Microgravity affects the physiology of organisms and the way they interact with each other and their environment. This may result in increased risk of infections.
- ▶ MuRGE is a flight experiment to study the effect of microgravity on the interaction between a plant (*Arabidopsis thaliana*) and a fungus (*Piriformospora indica*).



*Arabidopsis thaliana* (thale cress)



*Piriformospora indica*  
(green hyphae attached to root)




# MuRGE

- ▶ Why study plants and fungi in space?
  - Microorganisms (bacteria, archaea, fungi, viruses) are part of the spacecraft human environmental control and life support system. The microbes go where we go.
  - Plants can be a part of the life support system as well to provide fresh food, improve air quality (removing CO<sub>2</sub> and VOCs while providing O<sub>2</sub>), and recycle wastewater.
  - Living organisms (whether human, plant or microbe) are affected by living in reduced gravity. In some cases, this can result in an increased risk of infection and disease.
  - In normal gravity, the fungus *P. indica* increases the rate of plant growth and resistance to environmental stress.



# MuRGE


- ▶ Clinostats
  - Simulate the effects of reduced gravity on organisms by rotating them at low speed
- ▶ Controlled Environment Chamber (CEC)
  - 18:6hr Light:Dark, 400ppm CO<sub>2</sub>, 20°C, and 50% RH (relative humidity)





# MuRGE

- ▶ 36 agar plates (18 seeds/plate)
  - 12 horizontal, 12 vertical, 12 control
    - 3 At only
    - 3 Pi only
    - 3 At+Pi spores
    - 3 At+Pi extract




At and Pi extract (7x)– Day 4



# MuRGE

- ▶ Null hypothesis: *P. indica* will increase the growth of *A. thaliana* in both normal gravity and simulated microgravity.
- ▶ 2 experiments
  - Contaminated
  - No Pi growth
- ▶ Pi did not increase plant growth in normal gravity or simulated gravity





# Exposure to...

## Biolog

- ▶ Microbial species identification
  - 96-well plate
  - Yeasts, fungi, bacteria
  - Transmittance
  - Incubator– 30/37°C



## Other

- ▶ DNA extractions
- ▶ Acridine Orange (AO) Stains
- ▶ Microscopy
- ▶ Spectroscopy





# Work Related Activities

- ▶ Seminars
- ▶ VAB
- ▶ Intern socials





# Non-Work Related Activities




- ▶ Harry Potter Land
- ▶ Beach
- ▶ SAK comedy club






# Lessons Learned

- ▶ Be courteous and on time
  - ▶ Be flexible
  - ▶ Be prepared for anything
  - ▶ Don't be afraid to ask questions
- 



# Lasting Impact


- ▶ Lab skills
- ▶ Work experience
- ▶ Appreciation for NASA
- ▶ Understanding of workplace ethics
- ▶ Preparation for college





# Future Plans

- ▶ Georgia Institute of Technology
- ▶ Major: Biology
- ▶ Future job at NASA





# Acknowledgements

- ▶ Michael Roberts
  - ▶ Michele Birmele
  - ▶ Priscilla Moore
  - ▶ Amber Wade
  - ▶ Jim Gerard
  - ▶ Luke Roberson
  - ▶ NASA INSPIRE
- 



INSPIRE Pre-College Internship  
Kennedy Space Center  
Mutualism in a Reduced Gravity Environment (MuRGE)  
Patel, Karishma  
July 29, 2010

Reviewed by:

Michael S. Roberts, Ph.D.

Mail Code DYN-3

Dynamac Corporation

Engineering Directorate, Surface Systems Office (NE-S)

---



## **Abstract**

MuRGE (Mutualism in a Reduced Gravity Environment) is a NASA flight-research experiment to investigate the microgravity effects associated with cell-cell communication and beneficial microbe-host interactions using a plant-fungal model system. This investigation will use a clinostat, an instrument that slowly rotates the plants to negate the effects of gravitational pull on plant growth (gravitropism) and development, to simulate microgravity. I will be using the endophytic fungus *Piriformospora indica* (Pi) and the model plant species *Arabidopsis thaliana* (At). *P. indica* has been shown to colonize roots of various plant species, including *A. thaliana*, and to increase plant growth and resistance to stress. The fungus has the ability to grow from spores or in axenic cultures without the presence of a host. *P. indica* spores and *P. indica* extract will be used to inoculate *Arabidopsis* seeds germinated on a clinostat in order to determine if simulated microgravity affects the interaction between the fungus and its plant host.

## **Introduction**

Mutualism can be described as interactions between different organisms that are beneficial for both partners. One common effect of mutualism is increased growth of one or both organisms by increasing access to nutrients. Microgravity affects the physiology of organisms and the way they interact with one another and their environment. Plants and fungi go wherever humans go and are therefore a part of spacecraft environmental control and life support systems. It is necessary to study the effects of microgravity on microorganisms so we can be prepared for long duration missions in space. MuRGE will test the interaction between At and Pi in both normal gravity and microgravity environments. The null hypothesis is that *P. indica* will increase the growth of *A. thaliana* in both normal gravity and simulated microgravity.

## **Materials and Procedures**

This experiment will germinate seeds from *A. thaliana* (L.) Heynh wild-type ecotype Columbia on a solid agar media consisting of 2.2 g of Murashige and Skoog salts (Murashige and Skoog, 1962), 0.5 g of MES buffer, 2.5 g of sucrose, and 1 mL of 1,000x Gamborg vitamins (Sigma-Aldrich, St. Louis) per liter at pH of 5.75. Phytigel (Sigma) is added to a concentration of 0.8% (w/v) and autoclaved. After autoclaving, 5 mL of the media is aliquoted into a sterile 60 mm Petri plate (Millipore Microbiological Dishes, 47mm). Each Petri plate will be planted



with 18 sterilized *Arabidopsis* seeds per plate oriented in a single row across the plate on the solid surface. For *P. indica* spore treatments, *A. thaliana* seeds are coated with *P. indica* spores immediately prior to planting by immersing the seeds for 2 minutes in ~50  $\mu\text{L}$  of a spore suspension containing  $\sim 1 \times 10^4$  spores  $\text{mL}^{-1}$  dI  $\text{H}_2\text{O}$ . For *P. indica* extract treatments, *A. thaliana* seeds are coated with *P. indica* extract immediately prior to planting by immersing the seeds for 2 minutes in ~50  $\mu\text{L}$  of the extract.

Four slow-rotation clinostats rotating at 1 rpm, each holding 12 60-mm Petri plates arranged in 4 sets of 3, will be used during the experiment. Three Petri plates will contain only *Arabidopsis thaliana* (Treatment A), three plates will contain only *P. indica* (Treatment B), three will combine *A. thaliana* and *P. indica* spores (Treatment C), and a third set of plates will contain *A. thaliana* and an extract of *P. indica* (Treatment D). A third set of 12 60-mm Petri plates will be placed on the shelf as a no-centrifugal force control. Each Petri plate will be prepared in duplicate and one plate corresponding to each treatment will be covered in aluminum foil as an etiolated (dark) treatment. Two clinostats will be oriented vertically to rotate the Petri plates perpendicular to the Earth's surface (Simulated Microgravity) and the other clinostats will be oriented horizontally to rotate the Petri plates parallel to the Earth's surface (Normal Gravity Control). Petri plates loaded into the clinostats and the static controls will be placed in a Controlled Environment Chamber (CEC 15) under the following growth conditions: 18:6 Light:Dark,  $300 \mu\text{mol m}^{-2} \text{sec}^{-1}$ , 400 ppm  $\text{CO}_2$ ,  $20^\circ\text{C}$ , and 50% RH (relative humidity).

The growth of the plants will be monitored by measuring the length of roots and shoots of the *Arabidopsis* plants on the Petri plates and measuring the fresh-weight biomass. A single Petri plate will be harvested from each treatment at each of three four-day time intervals spanning twelve days. At each three-day interval, one plate from each treatment will be removed from the clinostat and the plants harvested for measurement under an Olympus SZX-12 zoom microscope.

#### **TREATMENTS (18 seeds per Petri plate; 3 Petri plates per treatment)**

A. <i>A. thaliana</i>	CTL ; HC ; VC	3 plates/treatment = 9 plates
B. <i>P. indica</i>	CTL ; HC ; VC	3 plates/treatment = 9 plates
C. <i>A. thaliana</i> and <i>P. indica</i> spores	CTL ; HC ; VC	3 plates/treatment = 9 plates
D. <i>A. thaliana</i> and <i>P. indica</i> extract	CTL ; HC ; VC	3 plates/treatment = 9 plates
		36 plates incubated in light
		36 plates incubated in dark



	<u>DOE4</u>	<u>DOE8</u>	<u>DOE12</u>
Control	CTL: A1/B1/C1/D1	CTL: A2/B2/C2/D2	CTL: A3/B3/C3/D3
Horizontal Clinostat HC:	A1/B1/C1/D1	HC: A2/B2/C2/D2	HC: A3/B3/C3/D3
Vertical Clinostat VC:	A1/B1/C1/D1	VC: A2/B2/C2/D2	VC: A3/B3/C3/D3

## Results

### 1<sup>st</sup> Experiment-Day 4

Plate			Orientation					Notes
	# of seeds	Germinated	Up	Down	Left	Right	Unknown	
A1 HC	18	14	6	4	2	2		
A1 N	18	17	15	2				
A1 VC	18	13	2	10	1			
B1 HC	18							16 spots of Pi (fungal contaminant?)
B1 N	18							Cannot distinguish individual spots, fungal growth along center
B1 VC	18							Cannot distinguish individual spots, fungal growth along entire line, one large growth on right
C1 HC	18	17	7	2	4	4		Spores and fungal growth not visible
C1 N	18	15	10	3	2			14 seem to have spores, 2 of those overtaken by Pi
C1 VC	18	16	5	5	3	3		11 visibly affected by Pi, others not visible
D1 HC	18	15	5	3	4	3	3	Plate had bubbles, unable to determine if affected by Pi
D1 N	18	13	11			2		Longest roots, all affected by Pi, 1 overtaken by Pi
D1 VC	18	16	3	10	1	2		Fungal growth visible across entire line, 1 overtaken



The data above was collected from the first set of plates, which contained eighteen seeds each. In most cases, more than three-fourths of the seeds germinated. In the neutral plates, a majority of the plants sprouted upwards. However, a few of the seeds that had recently germinated grew in other directions. The plates that were placed on the horizontal clinostat had similar results to one another. A majority of plants in each plate grew upwards. The plates on the vertical clinostat had a majority of seeds that grew downwards. The Pi growth was undistinguishable from the fungal contamination.

### 1<sup>st</sup> Experiment-Day 8

Plate			Orientation					
	# of seeds	Germinated	Up	Down	Left	Right	Unknown	Notes
A2 HC	18	17	14	3				
A2 N	18	17	17					
A2 VC	18	18	4	9	8	2		
B2 HC								10 visible colonies of contaminant (?)
B2 N								2 types of fungal contaminant
B2 VC								2 types of fungal contaminant
C2 HC	18	17	5	2	4		6	Pi/ contaminant has overtaken At- unable to observe roots/orientation
C2 N	18	16	6	5			5	Pi/ contaminant has overtaken At- unable to observe roots/orientation
C2 VC	16?	14	4	5	2	3		Plate contaminated, 5 overtaken
D2 HC	18	17						Plate dropped? All plants are loose and floating
D2 N	18	17	17					Bubbles around left side of plate, contaminate on right side
D2 VC	18	18	8	5	1		4	Many are shriveled, germinated at dif. times, newly germinated have been overtaken by contaminant

The data above was collected from the second set of plates from the first experiment. All of the plates, except one, showed to have eighteen seeds. A majority of those seeds germinated. In the *Aribidopsis* only and *Aribidopsis* + *P. indica* extract, all of the germinated seeds grew upwards.



However, the *Arabidopsis*+ *P. indica* spores had a randomly distributed amount of plants in all directions. The plates that were placed on the horizontal clinostat had a majority of plants that grew upwards. One plate had floating seeds, which could have resulted from the plate being dropped. The plates on the vertical clinostat had a randomly distributed amount of seeds growing in all directions. The direction and rate of Pi growth was unable to be quantified because of the presence of fungal contaminant(s).

## 2<sup>nd</sup> Experiment-Day 4

Plate			Orientation					
	# of seeds	Germinated	Up	Down	Left	Right	Unknown	Notes
A1 HC	18	18	9	4	4	1		
A1 N	18	15	15					
A1 VC	18	15	5	10				
B1 HC								1 distinguishable spot
B1 N								Pi not visible
B1 VC								Pi not visible
C1 HC	18	15	7	6		2		Pi not visible
C1 N	18	16	16					Root hairs visible
C1 VC	18	15	5	8	2			
D1 HC	18	14	7	4			3	A couple of roots are interwoven-unable to determine orientation
D1 N	18	18	18					Longest roots, root hairs somewhat visible
D1 VC	18	15	6	6		3		Pi not visible

The data above was collected from the first set of plates from the second experiment. In most cases, more than three-fourths of the seeds germinated. In the neutral plates, all of the seeds that germinated grew upwards. The plates that were placed on the horizontal clinostat had a majority of the plants grew upwards. The plates on the vertical clinostat had a majority of seeds that grew downwards. The Pi did not show any visible growth on these plates.



## **Discussion**


MuRGE investigated the microgravity effects associated with cell-cell communication and beneficial microbe-host interactions using a plant-fungal model system. This investigation used a clinostat to simulate the effects of microgravity. The plant *Arabidopsis thaliana* and the fungus *Piriformospora indica* were used to investigate these microbe-host interactions.

Two experiments were completed using the clinostats. While collecting data from the day four plates, a contaminant was spotted on the plates. The spores were the likely source of this. The inconclusive experiment resulted in starting a second experiment. Unfortunately, there was no visible spore growth in the first eight days of the experiment. Due to time restraints of my internship, I was unable to collect data for the day twelve plates of my second experiment. As a result of two inconclusive experiments, I was unable to determine the affects of Pi on At in both normal gravity and microgravity.


## **Acknowledgements**

I would like to thank my mentor, Michael Roberts, for all the time and support he gave me this summer and Michele Birmele for teaching me lab techniques and procedures. I would also like to thank all of the other staff at the Space Life Sciences Lab for making my summer experience an amazing one.











CEC 15



CAUTION  
CHECK CO2 LEVEL  
BEFORE  
ENTERING CHAMBER

CALIBRATE  
BEFORE USE

OUT  
13.0  
5.0  
SSR  
14.06

TC2 TEMP & RH BOARD  
ID (PLOC 20)=12

OUT  
12.50  
1.16  
SSR  
14.06

TC2 CO2 & LIGHT BOARD  
ID (PLOC 20)=13

SSR 1=LIGHT CAP 1 ENABLE  
SSR 2=LIGHT CAP 2 ENABLE

▲  
▼  
EGC  
TC2  
ENTER

▲  
▼  
EGC  
TC2  
ENTER

ON  
OFF

ALARM

RESET

OFF

RASTER CONTROL

OFF

LIGHTCAP #1

LIGHTCAP #2

Environmental  
Growth Chambers  
PO Box 390, Chagrin Falls, Ohio 44022  
USA



CEC 15



CAUTION  
HIGH VOLTAGE  
DO NOT ENTER  
ENTERING CHAMBER

CALIBRATE BEFORE USE

13.4  
5.0  
140.6

OUT

SSR RUN EDIT MOD

TC1 TEMP & RH BOARD  
ID: JPL0C 201112

12.45  
1.16  
140.6

OUT

SSR RUN EDIT MOD

TC2 TEMP & RH BOARD  
ID: JPL0C 201112

SSR T4 LIGHT CAP 1 ENABLE  
SSR T4 LIGHT CAP 2 ENABLE

EGC TC2 ENTER

Environmental Growth Chambers

800.321.8854  
440.247.5100  
440.247.5100 Fax  
www.egc.com

UNITCAP-1  
UNITCAP-2

GLOBAL CONTROLLED ENVIRONMENT CHAMBER # 36115

Temperature: 13.4°C  
Humidity: 50%  
Light: 140.6 μmol/m²/s

TC1 TEMP & RH BOARD  
ID: JPL0C 201112

TC2 TEMP & RH BOARD  
ID: JPL0C 201112

SSR T4 LIGHT CAP 1 ENABLE  
SSR T4 LIGHT CAP 2 ENABLE













































MEDIA 1 MEDIA 2



Age	Color	Example
1	Red	100% Red
2	Orange	100% Orange
3	Yellow	100% Yellow
4	Green	100% Green
5	Blue	100% Blue
6	Purple	100% Purple
7	Brown	100% Brown
8	Black	100% Black
9	White	100% White
10	Grey	100% Grey
11	Light Blue	100% Light Blue
12	Light Green	100% Light Green
13	Light Yellow	100% Light Yellow
14	Light Orange	100% Light Orange
15	Light Red	100% Light Red
16	Light Purple	100% Light Purple
17	Light Brown	100% Light Brown
18	Light Grey	100% Light Grey
19	Light White	100% Light White
20	Light Black	100% Light Black

